



Positioners



Our Set-Up technicians can process your Gnathological cases on virtually any articulator system, including the High Post Panadent™ and the Sam III™ instruments.



Our Positioner department has the skill to fabricate accurate and cost efficient "bench" Set-Ups and Positioners. We also specialize in partial and diagnostic Set-Ups.

The majority of the Positioner cases forwarded to AOA have the brackets reflected on the models. We accurately remove the brackets and fabricate the Set-Up and Positioner. This method allows the "Immediate" placement of the positioner upon clinical bracket removal, when teeth are most susceptible to guidance.



Pro-Visioner™ is made from casted Urethane to assure accuracy. Color and flavor choices are available.



Pro-Deucer™ and Pro-Former™ Urethane Positioners with enhanced freeway firmness in selected areas, posterior or throughout the appliance.



Flexiclear™ offers a wide variety of colors and flavors (see our Rx form) as well as great cost savings.



Accuracy, permanent flexibility and memory are the features of our Pro-Flex™ and Oralastic® silicone.



Im-Pak and Elast-Acryl™ are unique thermal sensitive materials that react to body temperature to produce a firm consistent force. These materials are ideal for extremely short trims – as little as 1/2 crown height.



Most of our adjuncts are available in all our materials, i.e., Claspings, Serrations, Hinge Axis Articulation, Air Holes and Sockette™ Reinforcements Liners.

COMMUNICATION CENTER

800-262-5221

INTERNATIONAL #262-886-1050

AOA Wisconsin

13931 Spring Street
Sturtevant, WI 53177
Fax 262-886-6879

AOA Connecticut

6 Niblick Road
Enfield, CT 06082
Fax 860-741-7655

AOA California

341 E. First Street
Calexico, CA 92231
Fax 760-357-9488

www.aolab.com

All of our Positioners and Mouthguards are custom trimmed to your exact specifications. The thickness and height can be tailored to produce the perfect blend of patient comfort and appliance efficacy. Custom Mouthguard information is obtainable by calling your Customer Service Representative.

What to send to AOA

Upper and lower models with suitably trimmed bases in either quality orthodontic plaster or stone and a completed AOA Rx form. The heels of the models should be trimmed into occlusion or a wax bite should be enclosed. If the brackets are to be left on, the archwires should be removed prior to taking the impressions.

AOA will return.

Set-Up models and Positioners with carrying case and patient pamphlet covering wear and care instructions. Your original models upon request for a small fee. Or, control models at no charge (simple untrimmed models without bases).