


# FROM ALLESEE ORTHODONTIC APPLIANCES

## Fixed Metal Appliances


AOA can fabricate almost any custom design of fixed appliance your office requires. We have chosen just a few of our most popular designs and modifications.



Kiddy Partial – one to four pontics can be provided. Just request the shade needed and a counter model if indicated.



Standard Lingual Arch – the anterior arch can be left ideal allowing fixed appliances to align the teeth. Adjustment "Omega" loops can be included upon request.



Habit Appliance – soldered crib at a vertical or slightly tipped angle to the patient's palate. Can be useful for either tongue or thumb habit.

Wire size and accessories are yours for the requesting...or call ahead with special questions or concerns. Soldered stops, springs, acrylic bite planes and pontics are all possible.


Quality plaster or stone models with seated bands are needed. AOA can supply bands if requested.

Additional information is available upon request to our Customer Service number:

**1-800-262-5221**

**Fax 1-414-886-6879**


[www.ormco.com/aoaapps.html](http://www.ormco.com/aoaapps.html)



Nance Holding Arch – also available with adjustment loops in the wire. The acrylic buttons can be requested in specific sizes by drawing the shape on the working model.



The Teflon Roll Bluegrass Appliance has been widely regarded as an effective tongue retraining aid.



Trans-Palatal Arch – soldered to bands or removable. Available in round or rectangular wire.



Fixed Bite Plane – this versatile appliance can be used for bite opening. Or, with an angle added to the acrylic, the appliance serves as an anterior repositioning fixed splint.

**sds AOA**  
SYBRON DENTAL SPECIALTIES

P.O. Box 725  
Sturtevant, WI 53177

AOA is a subsidiary of Ormco Corporation

Rev. B