



**Allesee Orthodontic Appliances**

**COMMUNICATION CENTER**

**1-800-262-5221**

**AOA Wisconsin**

13931 Spring Street  
Sturtevant, WI 53177  
Fax 262-886-6879

**AOA Connecticut**

6 Niblick Road  
Enfield, CT 06082  
Fax 860-741-7655

**AOA California**

341 E. First Street  
Calexico, CA 92231  
Fax 760-357-9488

INTERNATIONAL # 262-886-1050

**Flat Plane Occlusal Splint**

- Upper  Lower

Full coverage with even contact of opposing arch. Counter model suggested with wax bite. Used as bruxism or Pre Tx diagnostic splint.

**Retention Splint**

Requires upper and lower models and construction wax bite.

- Damon  
 Hard Pressure formed  
 Dual Hardness with Soft Liner  
 Elastic Silicone  
 Starnes Bite Orthotic  
 Moles Flexible Bionator

**Anterior Repositioning Splint**

- Upper  Lower

Full coverage upper appliance with anterior ramp designed to "capture" the lower anterior teeth and hold the mandible in a forward/anterior position. Lower appliance consist of posterior coverage indexed with upper lingual cusps connected with an anterior lingual band of acrylic. Forward upper and lower models. *Bite reflecting anterior position.*

**Bleaching Tray**

- Upper  Lower  
 Formed lined - .020 only  
 Vinyl no Bleach Reservoirs  
 Vinyl with Bleach Reservoirs

*Indicate Reservoirs on Diagram*

**Superior Repositioning Splint**

- Upper  Lower (Tanner splint)

Full coverage on upper arch with anterior ramp designed to provide protrusive and cuspid protection. Both upper and lower models with *centric relation bite required*. Bite should provide 1-2mm minimum anterior opening for acrylic durability.

Mounted models for adjustable articulators are recommended.

**MORA - Gelb**

Lower splint with posterior coverage connected with a heavy metal bar lingual of the incisors. Counter model is recommended with a wax bite reflecting desired occlusal acrylic thickness.


  **Check Box if using TAD's with Splint or Retainer**

**Draw TAD position on model**

**Anterior Open Bite Splint-AOB**

Upper posterior coverage with transpalatal bars. Buccal hooks are incorporated. Counter model advised. This appliance is used with implant screws (not available through AOA), which provides accelerated intrusion of upper posterior quadrants.

**Anterior Open Bite Retainer**



**SPLINT R<sub>x</sub>**

Dr. \_\_\_\_\_ Acct. # \_\_\_\_\_

Address \_\_\_\_\_

City St Zip \_\_\_\_\_

Patient \_\_\_\_\_


Tel \_\_\_\_\_ Fax # \_\_\_\_\_

E-mail \_\_\_\_\_

Shipped \_\_\_\_\_ Placement Date \_\_\_\_\_

**(PLACEMENT DATE SHOULD BE 1-2 DAYS BEFORE ACTUAL INSERTION DATE)**

PLEASE WRITE SPECIAL INSTRUCTIONS



**Materials**

- Splint Grade Acrylic  
 Pressure formed Biocryl base –  Add Acrylic  No Acrylic  
 Soft Vinyl appliance –  2mm  3mm  
 Dual Hardness – Soft inner surface with hard outer shell

**Color**  Clear Color \_\_\_\_\_

**Clasps**

Mark location on diagram

- Ball  Adams  
 Arrow  "C"

**Wax Construction Bite Provided**

- Use as is  
 Lab may modify if needed

**PLEASE SHIP EXTRA:**

- SHIPPING BOXES  
 PRE-PAID BAGS  
 PRESCRIPTION SHEETS

*Laboratory Use Only*

**SPL** \_\_\_\_\_

**SPL** \_\_\_\_\_

**ACRY** \_\_\_\_\_

**AAC** \_\_\_\_\_

**AAC** \_\_\_\_\_

**AAC** \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**SHIP** \_\_\_\_\_

**RECEIVING**

**OPEN** \_\_\_\_\_

**ID** \_\_\_\_\_

**DATE** \_\_\_\_\_

**ENTER** \_\_\_\_\_

**PULL** \_\_\_\_\_

**SHIPPING**

**PULL** \_\_\_\_\_

**LAYOUT** \_\_\_\_\_

**PACK** \_\_\_\_\_

**SHIP** \_\_\_\_\_

**CHECK** \_\_\_\_\_